

Detection of mineral nodules in marine fish mineralogenic cell lines

PREAMBULE. Mineralization of the extracellular matrix (ECM) of fish cell lines is previously induced by supplementing culture medium with 4 mM calcium chloride, 10 mM β -glycerophosphate and 50 mg/ml of ascorbic acid (Pombinho et al., 2004). Volumes and quantities given below are for a well of a 24-well plate.

PHOSPHATE-SPECIFIC STAINING (adapted from von Kossa, 1901)

1. Discard medium and wash cell culture 3 times with 1 ml of PBS.
2. Fix cells for 1 h at 4°C with 0.5 ml of 4% (v/v) formalin (in PBS).
3. Discard fixative and wash cell culture 3 times with 1 ml of MilliQ water.
4. Stain mineral nodules for 15 min at room temperature and under UV light with 0.5 ml of 5% (w/v) silver nitrate solution (in MilliQ water).
5. Discard stain and wash cell culture 3 times with 1 ml of MilliQ water.

Optional: Fix silver staining for 5 min at room temperature with 0.5 ml of 2.5% (w/v) sodium thiosulfate (in MilliQ water). Discard fixative and wash cells with 1 ml of MilliQ water.

6. Observe black-stained nodules under an inverted microscope equipped with phase contrast, then scan the plate at a resolution compatible with subsequent densitometry analysis (e.g. 600 dpi).

CALCIUM-SPECIFIC STAINING (ALIZARIN RED S / AR-S)

1. Discard medium and wash cell culture 3 times with 1 ml of PBS.
2. Fix cells for 1 h at 4°C with 0.5 ml of 4% (v/v) formalin (in PBS).
3. Discard fixative and wash cell culture 3 times with 1 ml of MilliQ water.
4. Stain mineral nodules for 15 min at room temperature with 0.5 ml of 40 nM AR-S solution (in MilliQ water).
5. Discard stain and wash cell culture 4 times with 1 ml of MilliQ water.
6. Observe coloration under an inverted microscope and document with micrographs and/or a plate scan.
7. Discard water and distain for 15 min at room temperature with 250 μ l of 10% (w/v) cetylpyridinium chloride (CPC) solution.
8. Collect extracts and measure absorbance at 550 nm (dilute if necessary). Determine stain concentration using AR-S standard curve (in CPC).

Vincent Laizé, Marta S. Rafael,
Sara Mira, M. Leonor Cancela
Centre of Marine Sciences (CCMAR),
University of Algarve, Faro, Portugal

Apparatus: An inverted microscope equipped with phase contrast, a UV light transilluminator, a spectrophotometer.

Solutions: phosphate-buffered saline solution (PBS: 137 mM NaCl, 2.7 mM KCl, 15.8 mM Na_2HPO_4 , 1.23 mM KH_2PO_4 ; adjust pH to 7.4); Formalin solution (4% in PBS; dilution from 36% commercial solution); AR-S solution (40 nM alizarin red S in MilliQ water; adjust pH to 4.2 with 1% NH_4OH). CPC solution (10% in 10 mM sodium phosphate pH 7.0).

Plasticware: 24-well cell culture dishes and serologic pipettes.

All chemicals were purchased from Sigma-Aldrich, unless otherwise stated.

Additional information:

Pombinho AR, Laizé V, Molha DM, Marques SMP, Cancela ML (2004) Development of two bone-derived cell lines from the marine teleost Sparus aurata; evidence for extracellular matrix mineralization and cell-type-specific expression of matrix Gla protein and osteocalcin. Cell Tissue Research 315: 393-406

von Kossa J (1901) Uber die im organismus kunstlich erzeugbaren verkalkungen. Beitr Path Anat 29:163-202