




Axenic *Ectocarpus* cultures

Axenisation of an *Ectocarpus* strain

1. Place antibiotic discs on Zobell medium with a flamed forceps (4 per dish 60 mm, 8 per 90 mm dish or 2 x 8 per 140 mm dish).
2. Place small fragments of *Ectocarpus* around and between the discs and seal the plate with parafilm.
3. Incubate at 18°C and 13°C and monitor the plates daily.
4. Take fragments of *Ectocarpus* from the areas where there is no development of bacteria with flamed forceps and put them in a 55 mm Petri dish containing natural seawater without Provasoli (to limit contamination).
5. Incubate the Petri dishes for 3-4 weeks at 13°C in a plastic box with a top or equivalent to limit contacts with Petri dishes containing non-axenic material.
6. Place a fragment on Zobell medium again, this time without the antibiotic discs. If no bacteria grow, the material can be cultured in natural seawater without Provasoli to produce biomass. If bacteria do appear, repeat the antibiotic treatment.

Plates with antibiotic discs can be kept for 1 month before carrying out the test. If large amount of biomass are needed, the material can be grown in a 10 litre bottle in PES for 1 month (without changing the medium, to limit the risk of contamination).



Müller DG, Gachon CMM, Küpper, FC

Fachbereich Biologie der Universität,
Germany

Institut de Biotechnologie des Plantes,
CNRS–Université Paris–Sud, France

Dunstaffnage Marine Laboratory, United
Kingdom


Additional information:

Müller, DG, Gachon, CMM, Küpper FC
(2008) Axenic clonal cultures of
filamentous brown algae: initiation and
maintenance. *Cah. Biol.Mar.* **49**: 59-65

Verifying that a strain is axenic

Working under a horizontal laminar flow hood, place some *Ectocarpus* material onto Zobell medium in a 60 mm Petri dish and incubate at 18°C (and, if possible, also at 13°C) for between 2 days and 4 weeks. Score by eye and under the microscope.

Bacterial contamination



The original material can also be observed under a microscope to look for the presence of bacteria.



Bacterial contamination

SOLUTIONS

ZOBELL MEDIUM (bacteriological solid medium)

Reagent	Quantity (for 500 ml)
Bactotryptone (Tryptone caséine triptone)	2.5 g
Yeast extract	0.5 g
Filtered natural sea water	400 ml
Distilled water	adjust the volume to 500 ml
agar (for solid media)	5 g

Pour into a 1 litre glass bottle and autoclave for 20 min at 120°C

Allow to cool slightly for 15-30 mn and then pour the culture medium into Petri dishes under the hood (if you don't pour the medium into the Petri dishes after 30 min, store the bottle in the oven at 60°C until use),

Allow the plates to cool with the lids open.

Seal with parafilm and store upside down at 13°C.

MATERIALS

Antibiotic discs

Antibiotic	Antibiotic/disc (μg) (Bio-Mérieux discs)	Colour (disc with initial)
Ampicillin	10	white
Chloramphenicol	30	white
Erythromycin	15	white
Kanamycin	30	white
Penicillin G	10	white
Rifampicin	30	orange
Tetracycline	30	yellow
Vancomycin	30	white

Axenisation of Ectocarpus cultures using the method published by Müller et al. (2008).